


TEETH WHITENING STRIPS

Best Teeth Whitening Products

PAIN FREE



SLS Free




No Animals Testing



100% Vegan



Peroxide Free



Fluoride Free

PRODUCT PARAMETER



Brand: Beautiful Smile

Product name: teeth whitening strips

Specification: 14 pouches/box, two week supply

Shelf life: 24 months

Usage: 1 time a day, 30mins per time

Target user: Food and drink, aging, smoking and other causes of tooth enamel exogenous yellow teeth

PRODUCT ADVANTAGES

- ✓ Mild formula
- ✓ Thin and invisible
- ✓ No pain no sore
- ✓ Non-slip design
- ✓ 360 degree whiten teeth
- ✓ 30 mins "teeth whitening mask"



PRODUCT INGREDIENTS



Menthol:
Fresh breath, cool and icy



Glycerin:
Moisturizing, moisturize, high activity, promoting alcoholization



Citric acid:
Gentle whitening, brighten teeth





Mild formula, completely safe. No side effect, no sensitivity. Do not damage the gums, protect tooth enamel, gentle brighten teeth.

SUGGEST 1+1 WHITENING TREATMENT 28 DAYS

1+1 whitening treatment:

1 box basic whitening + 1 box consolidate whitening = 1 treatment

After 28 days, use 2-3 times per week to last the whitening effect.




- ✓ 7 days: give you a surprise smile in 7 days
- ✓ 14 days: show your dazzling white teeth in 14 days
- ✓ 28 days: you can control any color lipstick in 28 days

14-DAYS WHITENING TEST


Brightness improvement degree




MAKE UP YOUR TEETH AT IMPORTANT MOMENTS



Interview



Wedding



Dating

Small and portable, Does not occupy bag space
Lipstick CP, Red lips and white teeth are more attractive



EASY OPERATION NO DIFFICULTY FOR NOVICE





Step 1: Peel

Open strips packet, tear off the strips from the backing liner

Step 2: Apply

Apply gel side against your teeth, lining up with edge of gum. Apply gentle pressure to mould the strips securely to the teet



Step 3: Reveal

Peel off strips after 20-30 mins, rinase mouth with water.

OPTIONAL FORMULA



6% HP



Non Peroxide



Activated Charcoal



PAP

CUSTOMIZATION



Custom pouch

Custom box